

Communication Solutions

From simple to complex, speech generating devices (SGDs) for augmentative and alternative communication (AAC) assist persons with disabilities at all levels of language development starting with simple, single message devices. Users can easily transition to basic choice making, sequential communication, and even more complex communication opportunities. AbleNet's complete line of SGDs is there for you at each stage.

When choosing a SGD, a few questions to consider:

1. How many messages will be required?
2. Where will the device be used — on a table, on a wheelchair, or in hand?
3. What access methods will be used, direct or external switches?
4. Does this device need mounting?
5. Will this device be used by an individual or in a group setting?

Communication Solutions												
MESSAGES / LEVELS		1 / 0	1 / 0	Unlimited sequential / 3	2 / 3	4 / 3	Unlimited sequential / 3	4 / 2	7, 12, 23 / 5	1, 2, 4 or 8 / 8	Varies by App	Varies by App
RECORDING TIME		10 sec	2 min	4 min	4 min	6 min	4 min	80 sec	6 / 12 / 20 min	16 min	Varies by App	Varies by App
SWITCH INPUT		No	1	1	1	4	1	No	No	8	Includes Blue2™	Includes Blue2™
TOY / APPLIANCE OUTPUT / ECU		No	1	1	1	4	1	No	No	2	No	No
ACTIVATION SURFACE (in/cm)		1.8 / 4.6	B: 5 / 12.7 L: 2.5 / 6.4	B: 5 / 12.7 L: 2.5 / 6.4	2.5 / 6.4	2.5 / 6.4	2.5 / 6.4	0.5 x 0.7 / 1.3 x 1.8	Various sizes 1.2 x 1.2 / 3 x 3 to 3.4 x 2.7 / 8.6 x 6.9	Various sizes 1.8 to 7.8 / 4.5 to 19.8	Varies by App	Varies by App
ACTIVATION FORCE (oz/g)		2 / 57	B: 5 / 142 L: 3 / 85	B: 5 / 142 L: 3 / 85	5 / 142	5 / 142	3 / 85	5 / 142	3 / 85	0 / 0	0 / 0	0 / 0
WEIGHT WITH BATTERY (oz/g) BATTERY TYPE		1.5 / 43 Rechargeable	B: 15.5 / 439 L: 8.5 / 241 9-volt	B: 15.5 / 439 L: 8.5 / 241 9-volt	12.5 / 354 9-volt	23.2 / 658 4 AA	8.5 / 241 9-volt	0.12 / 3.4 Rechargeable	26.4 / 748 4 AA	33.6 / 953 4 AA	37 / 1,048 Lithium ion	27 / 765 Lithium ion
ACCESSORIES AVAILABLE		N/A	Travel Strap, Mounting Plate, Interchangeable Switch Tops	Travel Strap, Mounting Plate, Interchangeable Switch Tops	Travel Strap, Mounting Plate	Mounting Plate	Travel Strap, Mounting Plate, Interchangeable Switch Tops	Works with alternative wrist straps	N/A	Travel Strap, Mounting Plate, Adjustable Stand	Keyguard, Mounting, Shoulder Strap	Keyguard, Mounting, Shoulder Strap
MOUNTING CONNECTION		Internal Magnet	80000015	80000015	80000015	80000015	80000015	N/A	N/A	80000015	Latitude QTF	Latitude QTF
HCPC CODE		E2500	E2500	E2500	E2500	E2500	E2500	E2500	E2500, E2502, E2502	E2502	E2510	E2510
ITEM #		10000032	B: 10002100 L: 10002400	B: 10002300 L: 10002600	10003300	10000028	10003600	10002000	10003501 10003502 Not Available	10002800	70000117	70000119
COST		\$199	\$135	\$195	\$190	\$250	\$195	\$145	\$185, \$225, N/A	\$410	Call for pricing	Call for pricing