

RE: Ensuring the QuickTalker Freestyle can ONLY be used as a communication device

To Whom It May Concern:

The QuickTalker Freestyle iPad-based device is provided to assist with speech as a medical device for individuals who are non-verbal. The device and the associated application software allow an individual to press a selection on the screen and the result will be a word or phrase spoken for the individual. The device gives the user the opportunity to communicate his or her needs independently.

Although the device is iPad-based, it has been technically configured to ONLY allow access to the communication application feature. All other typical or common features found on an iPad have been disabled making the QuickTalker Freestyle a dedicated medical device for speech use ONLY. Therefore, functionality such as camera use, video recording, and sharing are unavailable to the user in the current state of the configuration. Because these features have been disabled, the device should not violate any policies that you have established for any of your customers regarding intellectual property at the venue.

To verify the above features have been disabled on the device, we encourage you to test the following:

- Taking a picture via the Camera app
- Taking a video via the Camera app
- Tapping the share icon available in most apps

If you have additional questions, please reach out to the ableCARE Product Success Team using the following options:

- Email: ableCARE@ablenetinc.com
- Phone: (651) 294-3101
- Via the ableCARE app on your QuickTalker Freestyle device

Thank you,



Mike Knott

Manager, Customer Relations

AbleNet, Inc.
2625 Patton Road, Roseville, MN 55113
Direct (651) 294-2225
Email mknott@ablenetinc.com
Web www.ablenetinc.com / www.quicktalkerfreestyle.com

Download the ableCARE app for access to on-demand support and resources to ensure best-outcomes with your AbleNet assistive technology and QuickTalker Freestyle speech device. Available on the [App Store](#) and [Google Play](#)